

بِسْمِ اللَّهِ الرَّحْمَنِ الرَّحِيمِ
اللَّهُمَّ زِدْنِي عِلْمًا

Tissue Culture Lab & Techniques

S	Date	Topic
Lab(1)		Tissue Culture Lab and Equipment (Laminar flow hood- Liquid nitrogen container-Autoclave-Centrifuge-Water bath-Ice maker- Inverted microscope-pipettes-flasks- How to use the pipettes (blue and yellow tips))
Lab(2)		Preparation of media (use of magnetic stirrer and pH meter) Sterilization (Filtration & Autoclaving)
Lab(3)		Thawing and growing Primary culture : Liver (adherent) and spleen (suspension) cells Growing cells on cover slips
Lab(4)		Cell count (Haemocytometer) Primary culture: monocyte cell (Method of Extraction & Culture)
Lab(5)		Freezing Fixing and staining cells grown on cover slips
Lab(6)		Photographing cells with digital camera and printing
Lab(7)		Presentation
Lab(8)		Exam

LAB (1)

Tissue Culture Lab and Equipment

Essential Equipment :

- Laminar Flow Hood (Safety Cabinet) :
Gives best sterile protection to the culture.
- CO₂ Incubator :
A device with uniform temperature can be maintained, used in growing cultures.
- Autoclave :
Sterilization under high temperature and pressure.
Solid materials sterilize in 121°C/15 min and liquid in 121°C/20 min
- Inverted microscope :
Quantitative analysis of cultured cells for seeing the sample in three dimensional.
- Magnetic Stirrer :
Rapid stirring action for dissolving chemicals is available with any stirrer for disaggregating or suspension culture.
- pH meter :
Adjusting the pH.
- Liquid Nitrogen Container :
To save the sample as it is, for long time in -198°C
- Sensitive Balance :
To weight the substances.
- Vacuum Air Pump :
A sterilizing by filtration (fluid pass throw filter membrane to remove all type of microorganism in the culture medium.
- Centrifuge :
There are two sizes
Microfuge which used with Eppendrof tubes.
Macrofuge which used with 3,5,10,15,50ml tubes.
- Refrigerator.
- Freezer(-20°,-80°)C.

- **Plastic Well Plate:** To seeds the cells.
- **Cry vial Tubes:** To freezes the cell in liquid nitrogen container.
- **Hot Plate.**
- **Homogenizer.**
- **Water Bath.**
- **Milex filter.**
- **Stir Cup filter.**
- **Shaker.**
- **Ice maker.**
- **Pipettes :**

There are different types:

Multi pipette – single pipette – simple pipette – automatic pipette – micro pipette – disposable (Pester) pipette.

Tips:

Single pipette:

Crystal tips: 0.5 -10 μ l

Yellow tips: 10-200 μ l

Blue tips: 200-1000 μ l

White tips: 1- 5 ml

- **Flask :**

To seed the cells.

There are two sizes: 75 cm³ & 25 cm³

When we put the flasks inside the incubator we must loosen the caps to allow sufficient gas exchange.

How to use pipettes :

Single pipette –automatic pipette – disposable Pester pipette.

Preparation of solutions with different concentrations

Solutions & dilutions :

$$C_1V_1 = C_2V_2$$

C_1 = starting concentration. V_1 = starting volume.

C_2 = final concentration. V_2 = final volume.

Want × Final volume

Have

w/w (weight / weight) solutions :

1% solution is 1g solute in 100g final weight of the solutions.

1ml H₂O weight 1g

Mix weight of solvent with volume of H₂O

- 10% NaCl = 10g NaCl in 90ml H₂O

w/v (weight/ volume) solutions :

1% w in v means:

1g solute in a final volume of 100 ml.

Examples :

*Prepare 50 ml of 2% NaCl.

2 g _____ → 100 ml H₂O
? _____ → 50 ml

$$\frac{50 \times 2}{100} = 1 \text{ g}$$

*Prepare 25 ml of 2% NaCl .

2 g _____ → 100 ml H₂O
? _____ → 25 ml

$$\frac{25 \times 2}{100} = 0.5 \text{ g}$$

*Prepare 200 ml of 0.85% NaCl.

0.85 g _____ → 100ml H₂O
? _____ → 200ml

$$\frac{0.85 \times 200}{100} = 1.7 \text{ g}$$

***Prepare 50 ml of 0.3% NaCl.**

0.3g _____ → 100 ml H₂O
? _____ → 50 ml

$$\frac{0.3 \times 50}{100} = 0.15 \text{ g}$$

v/v (volume/ volume) solutions :

Examples:

***Prepare 100 ml of 10% acetic acid.**

10ml acetic acid + 90ml H₂O

***Prepare 20ml of 0.2% acetic acid.**

0.2ml _____ → 100ml H₂O
? _____ → 20ml

$$\frac{0.2 \times 20}{100} = 0.04 \text{ ml}$$

***Prepare 50ml of 0.5% acetic acid.**

0.5ml _____ → 100ml H₂O
? _____ → 50ml

$$\frac{0.5 \times 50}{100} = 0.25 \text{ ml}$$

*** Prepare 70% alcohol.**

70ml alcohol + 30ml H₂O

*** Prepare 400ml (1x) from (10x) stock PBS.**

$$\frac{1}{10} \times 400 = 40 \text{ ml}$$

$$400 - 40 = 360 \text{ ml H}_2\text{O}$$

Molarity :

To prepare a solution in Moles/Liter = M

The molecular weight (mw) must be known.

E.g. mw of NaOH =40

40g in 1L (1000ml) =1M

Mw × concentration required = amount to weight/liter

$$40 \times 0.5 = 20\text{g in 1L}$$

Volume required = 500ml

Amount/liter × volume required = amount to weight

Examples:

*prepare 100ml of 0.5 M NaCL. (Mw= 58.88)

$$0.5 \times 58.88 = 29.44$$

58.88 In L in 1M

29.44 In L in 0.5

? _____ 100ml

$$\frac{100 \times 29.44}{1000} = 2.944 \text{ g}$$

*Prepare 200ml of 70mM sucrose. (Mw = 342.3)

$$70\text{mM} = \frac{70}{1000} = 0.07 \text{ M}$$

$$200\text{ml} = \frac{200}{1000} = 0.2 \text{ L}$$

$$342.3 \times 0.07 = 23.961$$

342.3 In 1L in 1M

23.961 In 1L = 0.07 M

? _____ 0.2

$$23.961 \times 0.2 = 4.7922$$

*Prepare 100ml of 50mM sucrose. (Mw = 342.3)

$$50\text{mM} = \frac{50}{1000} = 0.05 \text{ M}$$

$$100\text{ml} = \frac{100}{1000} = 0.1 \text{ L}$$

$$342.3 \times 0.05 = 17.115$$

342.3 In 1L in 1M

17.115 In 1L = 0.05 M

? _____ 0.1

$$17.115 \times 0.1 = 1.7115$$

Dilutions:

Prepare 400 ml (1 x) solutions from (10x) stock PBS.

$$C_1V_1 = C_2V_2$$

Examples:

❖ How much ml of 12 M HCl is required to prepare 250 ml of 2 M HCl solution.

$$\frac{2 \text{ M} \times 250}{12 \text{ M}} = 41.6 \text{ ml}$$

$$250 - 41.6 = 208.4 \text{ H}_2\text{O}$$

❖ What volume of 10 M acetic acid is required to prepare 1L of 0.50 M acetic acid?

$$\frac{0.050}{10} \times 1000 = 50 \text{ ml}$$

$$1000 - 50 = 950 \text{ ml H}_2\text{O}$$

❖ Dilute 500 ml 325 Mm sucrose to final concentration of 200 mM .

$$325 \text{ mM} = \frac{325}{1000} = 0.325 \text{ M}$$

$$500 \text{ ml} = 500 = 0.5 \text{ L}$$

$$C_1V_1 = C_2V_2$$

$$0.325 \times 0.5 = 0.2 \times V_2$$

$$\frac{0.325 \times 0.5}{0.2} = 0.812 \text{ L}$$

$$0.2$$

LAB (2)

Preparation of media&Stains

1) Preparation of RPMI media (1L):

- 1- Add 700ml distilled water to the ready measured powder media (15.8g or as it required in the container).
- 2- Add 20ml 10% sodium bicarbonate (NaHCO_3).
Prepare 100ml of 10% NaHCO_3 (10g NaHCO_3 +90ml dD.W)
Filter (NaHCO_3) (with milex filter 0.45 μm membrane)
- 3- Add 1% (10ml) L-Glutamine (200mM)
- 4- Add 1% (10ml) antibiotic
(penicillin10.000IU/streptomycin10.000 μg \ml).
- 5- Add 160 ml distilled water.
- 6- Adjust the pH 7.2 - 7.4.
If pH is high: add HCl (acid).
If pH is low: add NaOH (base).
- 7- Filtering through Millipore filter with 0, 22 μ pore size using vacuum air pump.
- 8- Add 100ml FCS (Fetal Calf Serum).

Preparation of stains:

A) Giemsa Stain (10%):

- 1- Add 10ml Giemsa Stain to 90ml distilled water.
(Filter with filter paper every time before using).

B) Trypan Blue (0.4%):

- 1- Add 0.4g Trypan Blue
- 2- 100ml Saline (0.9g NaCl in 1L D.W)
- 3-Filter (0.45 μm membrane).

LAB (3)

Growing Cell Line :

<u>Cell Line</u>	<u>Cell Type</u>	<u>Medium</u>
MCF-12	Adherent	Universal Media
MCF-7	Adherent	RPMI
MDA	Adherent	RPMI
HSFN-1	Adherent	DMEM-F12
THP-1	Suspension	RPMI

Thawing Cry preserved Cells:

- 1- Wear cry vial gloves and face shield and retrieve the required (check the index) ampoule (cry vial tube) from the liquid nitrogen container.**
- 2- Defrost the vial in a 37°C water bath with constant, moderate agitation.**
- 3- When thawed, swab the vial with 70% alcohol.**
- 4- Working in Sterile area (hood), open the vial and transfer the contents using sterile Pasteur pipette in to 15ml centrifuge tube(V shape).**
- 5- Add 10ml pre-warmed (37°C) culture medium gradually *.**
- 6- Centrifuge the cells 1500 rpm for 5min.**
- 7- Discard the supernatant and gently resuspend the cell pellet in 15ml culture medium.**
- 8- Transfer the cells in a 25cm³ culture flask.**
- 9- Observe the cells and place them in the CO₂ incubator.**

***this is important with Dimethyl Sulfoxide, where sudden dilution can cause severe osmotic damage and reduce cell survival by half.**

Primary Culture :

- Growing Liver Cells (adherent) .**
- Growing Spleen Cells (suspension) .**

Initiation of a Primary Culture :

- 1- Dissect the animal.
- 2- Take the liver and the spleen and wash with sterile normal saline.
- 3- Homogenize the organs with 10ml media (insert the homogenizer in ice).
- 4- Transfer the homogenized organs into 15ml centrifuge tubes and Close tightly.
- 5- Centrifuge 500 rpm for 5min.
- 6- Transfer the supernatant into a 15ml centrifuge tube.
- 7- Centrifuge 2500 rpm for 15 min.
- 8- Discard the supernatant and add 10ml media to the pellet shake well.
- 9- Take 1ml of the tissue and add to 9ml media.
- 10- Transfer to a labeled 25cm³ flask.
- 11- shake the flask gently and put in CO₂ incubator 37°C.

Notice :

The liver cells grow as an adherent monolayer.
The spleen cells grow as suspension.

Staining technique:

- Giemsa Stain.
- 1- Dissect a rat (or mice) and cut out a small piece of the liver and spleen.
 - 2- With a forceps take a small piece of both organs.
 - 3- Place a small amount of the specimens (by pressing down) on clean microscope slides.
 - 4- After the specimen in the slide dry fix it in 100% methanol for 1min.
 - 5- Stain with Giemsa stain for 25 min.
 - 6- Wash the slides with running tap water.
 - 7- Dry the slides with filter paper.
 - 8- View with the light microscope.

Sub culturing Adherent Cells:

- 1- Discard media from the culture flask.**
- 2- Rinse cells with PBS (10ml/75cm³- 3ml /25cm³) by gently rocking the flask back and forth, then remove the PBS and discard.**
- 3- Add trypsin to the flask (3ml/75cm³-1ml/25cm³) and rock the flask to ensure that the entire monolayer is covered with trypsin solution.**
- 4- incubate the cells 3-5 min until the cells detach (avoid leaving cells exposed to trypsin longer than necessary and to force cells to detach as this may result in clumping) .**
- 5- Add culture media (pipette the cells up and down) then transfer in a 15ml centrifuge tube.**
- 6- Centrifuge the cells 1500 rpm for 5 min *.**
- 7- Discard the supernatant and resuspend the cell pellet in fresh culture media, pipette the cells up and down.**
- 8- Keep the cells in ice.**
- 9- Count the cells with chamber (hemocytometer) then dilute to the appropriate concentration for seeding.**
- 10- Add the appropriate volume of cells suspension to a new flask containing media.**
- 11- Place flask in incubator with loosens caps to allow gas exchange.**

***(step 6 can be avoided as the serum in the medium will inactivate the remaining trypsin and prevent cell damage).**

Sub culturing Suspension Cells:

- 1- Transfer the cell suspension in a centrifuge tube.**
- 2- Centrifuge the cells 1500 rpm for 5min.**
- 3- Resuspend the cell pellet in culture medium.**
- 4- Count the cell by hemocytometer and dilute to the appropriate concentration for seeding.**
- 5- Add the appropriate volume of cell suspension to a new flask containing medium.**
- 6- Place flask (standing) in incubator with loosens caps to allow gas exchange.**

LAB (4)

Hemocytometer Counting Chamber:

- 1- Clean the surface of the hemocytometer slide and cover slip with 70% alcohol.
- 2- Place the cover slip over the hemocytometer counting chamber.
- 3- Shake the cell suspension well and mix equal volume of cell suspension and Trypan Blue stain (1:1).
- 4- Transfer 20µl cell suspension to the edge of the chamber and allow it to be drawn into the chamber by capillary action.
- 5- Place the hemocytometer on the microscope stage.
- 6- Count cells in the four marginal chambers using 40 X objective.
- 7- to calculate the number of cells :

$$c = n/v$$

c = cell concentration (cells/ml).

n = average number of cells counted.

v = volume counted = 10^4

thus : $c = n \times 10^4$

Concentration of cells/ml = average number of cells x DF* X 10^4

* DF: dilution factor

- 8- after calculating the cell concentration for seeding :

Concentration of cells required x final volume

Concentration of cell /ml

- 9- Add the appropriate volume cell concentration required to a new flask containing medium (final volume).
- 10-place flask in incubator with loosens caps to allow gas exchange.

LAB (5)

Preparation of the freezing media:

RPMI 1640 media	80 ml.
Serum	10 ml.
Dimethyl Sulfoxide (DMSO)	10 ml.

Freezing:

- 1- Count the cells.
- 2- Centrifuge the cells 1500 rpm for 5 min.
- 3- Discard the supernatant and add the required amount of freezing media to the pellet.
 - * Amount of freezing media needed to be added to the cells:
$$\frac{\text{Cell count obtained/ml}}{\text{Cell count required/ml}}$$

(The freezing media requires 10^6 /ml cells)
- 4- Distribute the cells in cry vial labeled tubes.
- 5- Place the tubes in freezing box (Iso propanol) -1°C per mints.
- 6- Place the tubes in the freezer (-80°C) for 3-4 days.
- 7- Transfer the tubes from the freezer to the liquid nitrogen container.

Growing Cells In 24 Well Plate

- 1- Subculture the cells.
- 2- Determine cell concentration.
- 3- Insert sterilized cover slips (13mm) in a 24 well plate using sterilized forceps and needle.
- 4- Calculate the seeding density required for each well :
$$\frac{\text{Cell concentration required} \times \text{Final volume}}{\text{Cell count/ml}}$$
- 5- Place 1ml of the cell concentration required in each well (label the plate)
- 6- Shake the plate gently, cover and place in CO₂ incubator 37°C.

Fixing & Staining the Cells with Giemsa - H& E stain

A) Fixing:

- 1- Remove media from the 24 well plate.
- 2- Fix the cells on the cover slip with 10% formalin (or methanol) for 10min.
- 3- Wash with PBS solution.
(Dissolve one tablet of PBS in 100ml sterile D.W)

B) Staining (Giemsa):

- 1- Stain the cells with Giemsa for 20 min.
- 2- Wash in tap water.
- 3- Remove the cover slips from the plate (wells) .
- 4- Prepare clean labeled slides.
- 5- Mount with DPX media.

C) Staining with (H&E):

- 1- Stain with Haematoxylin for 7-10 min.
- 2- Wash in tap water and rinse in distilled water.
- 3- Stain with Eosin for 1min.
- 4- Wash in tap water and rinse with distilled water.
- 5- Remove covers lips from wells.
- 6- Prepare clean labeled slides.
- 7- Mount with DPX.

LAB (6)

Photographing cells with digital camera and printing


Inverted Microscope




Compound Light Microscope



Image Management System



Archiving



Movie Recorder



*Movie Recorder/
Movie Player Module*



Image 1: Object in focus



Image 2: Background in focus



Image 3: Object and Background in focus

Measurement Module

